





SMARTCHEM 600

THE NEW GENERATION DISCRETE ANALYZER

SMARTCHEM

AMS Alliance

SMARTCHEM 600

Fully automated, direct read, discrete Analyzer to maximize productivity For routine and special chemistries, UP to **600 tests/hour**



AFTER HAVING SUCCESSFULLY INTRODUCED A COMPLETE RANGE OF DISCRETE ANALYZERS IN LABORATORIES, THE SMARTCHEM RANGE (OVER 7000 DISCRETE SYSTEMS INSTALLED WORLDWIDE), WE HAVE DEVELOPED THE SMARTCHEM 600 WHICH PROVIDES YOU:

MORE PRODUCTIVITY

LOWER RUNNING COSTS

BETTER RESULTS

The Smartchem 600 combines the latest technology in discrete analysis with an innovative and unique design that guarantees accurate and reproducible results.

The Smartchem range is the first discrete analyzers to utilize washable cuvettes with an integrated wash and control station, which guarantee a lower running cost and reduced risk of contamination. In addition to automated analysis the Smartchem 600 offers automated system Quality Control.

The Smartchem range is also the first discrete analyzers to utilize a patented, dedicated sample prep module for the determination of NOx by Cadmium reduction. Designed for multi-tasking and high throughput rates, the NOx module also provides automated cadmium regeneration. The instrument provides also a pH, Conductivity and Redox module.

Fully automated, the Smartchem 600 is available with new user-friendly software designed for touch screen computer with the possibility of wireless control.

Ease of use combined with minimal operating condition requirements make the Smartchem 600 the right choice for every cost-efficiency laboratory.

MAIN ADVANTAGES

- 600 tests per Hour
- 200 samples capacity Automatic startup and
- Automatic startup and calibration routine
- Unlimited number of unattended tests
- Barcode reading
- No disposable cuvette
- Integrated wash & control station
- Caped tubes
- New software with touch screen
- Bi-directional LIMS ready

MORE SIMPLICITY FOR BETTER RESULTS

The Smartchem 600 utilizes direct read photometry utilizing a true discrete, direct read measurement system.

- Fully automated operation, just load the Samples and start
- 6 mm or 10 mm Optical Path for increase accuracy
- Complete range of reagents and methods

MORE PRODUCTIVITY


- 200 Samples capacity (2 carousels with 10 racks)
- Continuous loading
- STAT sample
- Barcode reading for samples & reagents
- Repeat & Reflex tests
- Up to 100 reagents on board

LOWER RUNNING COST

- Limited space occupancy
- Use of washable cuvettes
- Reduced hazardous waste management cost
- Reduced volume of waste
- Low energy & Water consumption
- Green technology

ISE module for pH, conductivity & redox analysis










3 arms, one for samples, two for reagents



Long-life cuvettes with integrated Cuvette Quality Assurance System



Sample capacity 200, with continuous loading



Enabling use of tablet to monitor your analyzer

CUVETTE QUALITY ASSURANCE SYSTEM: LESS RISK OF CROSS CONTAMINATION



- Automated cleaning eliminates carryover in the cuvette
- We clean, rinse, dry and verify optically each cuvette before and after each single test
- Provides an infinite number of reaction vessels for analysis
- Economical: long lasting cuvettes

CAPPED TUBES



With the new capped tubes, you can run any type of samples without any problems.

- You can load low and high concentrated samples on the same rack without any risk
- The caps will prevent evaporation and contamination.

OPTIONAL MODULES

Multiple techniques are available to answer all environmental application and include electrochemical parameters:

- ISE module for pH, conductivity and redox analysis
- NO3 module with new reductor

NEW USER-FRIENDLY SOFTWARE WITH INTERFACE DESIGNED FOR TOUCH SCREEN COMPUTER

The Smartchem 600 use a new user-friendly operating software with intuitive access concept and touch screen approach. Extremely flexible, it is designed for either keyboard or touchscreen. It can be used with remote assistance service.

NEW SOFTWARE MAIN FEATURES

- Automatic start-up
- Automatic shutdown
- Barcode for reagents
- Barcode for samples
- Addition of new sample or method to a sample during the run
- Reanalyse key
- Reflex testing
- Accurate estimation of completion of analysis
- Calibration points recalculation
- Last accepted calibration overlay
- Pre-dilution
- Standard curves
- Ready to read 2D barcoding in reagents kits
- Computer controlled micro pipettor
- EPA Nilac compliant

MAIN PARAMETERS

Surface, ground, drinking, waste and sea water

- Acidity
- Alkalinity*
- Aluminum .
- Ammonia*
- Chemical Oxygen Demand (COD)* Chloride*
- Chlorine
- Wine
- Acetic acid •
- Ammonia
- Ascorbic Acid
- Cupper
- Ethanol (TVA)
- Free SO2

Soils & Plants

- Ammonia
- Boron
- Calcium
- Chloride

- Chromium*
- Cyanide* •
- Fluoride* •
- Hardness*
- Manganese
 - Nitrate*
- Nitrite*
- - Total Kjeldahl Nitrogen* Total Kjeldahl

Phenol*

Silicate*

Sulfate*

Phosphate*

Phosphorous*

рΗ

•

•


•


- Sorbic Acid •
- Sucrose
- Tartric acid
- Total acidity Total SO2

- - Magnesium

 - Nitrite
- Phosphate
- Total Nitrogen
- Total Phosphorous


- **Feeds & Fertilizers**
- Ammonia
- Calcium
- Chloride
- Phosphate
- Protein
- Total Nitrogen
- **Total Phosphorous**





0 _ 4		Straight and Torons		Renard B Serley		Quality Control			
the local and	~	colour and mean	Anger .	APOPT A SPAN	Canada	drawly research		INFO MALE CONTRACTOR	
National And		Result under trearily						19-29 - 504 - 0-20 20x 5-504	
. Mitch		P _ 27		12	10	18		18-495-194-0-1079115-1 water	
Artist	Refer	Alteriations 🖌						NEDA - 685 1580 NO27 mpt	_
Pulles								NEEL Correl - ANS \$142 ME22 mail	
Peterder	Reyn			17	10			1402Y - 1402 S ppm N	
Serging ode:	Refer	Abelet enr 🖌		2	1	1	1	NOW - NOCI por N	
B Chells	>			and high (2 mg/L)				https://dilace/t-leadow	
B faulter	3 500	1 J.				10		PDH - PD4 2-13 april 7	
🛃 Aulter		Martan ana		1	1	1	1	NORM-19040-63 ppn 7	
Ø Heistereret	>			and our leasily	-	-	-		
Options	>			T.					
in my	,	Delatana d		*	1		-1		
0805	>	And maximum 4					-		
3									
Xore								+ -	Crel

- - 1 Up to 100 reagents
 - 2 Multipoint calibration, Last
 - calibration, Compare & overwrite
 - Rerun & Reflex tests, Automatic 3 spiking recovery calculation
 - 4 200 samples at a glance
 - 5 touch screen computer



OD 420, OD 520, 0D 620 PANOPA

Malic acid

Gluconic acid Glucose

READY TO USE REAGENTS

- Developed specifically for the SMARTCHEM analyzers
- A full range of reagents for Environmental and Wine analysis
- Complete range adapted to all parameters
- Save time
- Save money
- Simple
- Ready to use



We have created Selection, a full range of enzymatic and colorimetric reagents dedicated to discrete analyzers Smartchem. Developed, tested and validated by our application laboratory, they will guarantee a perfect stability and reproducibility of your analyses. Very easy to use, they are supplied liquid, in specific bottles directly positionable on Smartchem, and automatically detected by the bar code reader (in option). Thus, this ready solution

greatly simplifies the preparation of your analysis, facilitates the use and storage of reagents. It also secures yours analyzes with reagents adapted to your application and your instrument.

REAGENTS FOR ENVIRONMENTAL ANALYSIS	REAGENTS FOR WINE ANALYSIS
Selection	Selection
Alkalinity Ammonia Chromium VI Chlorides Cyanides Nitrate Nitrate Nitrite Phosphates Silicate Sulfates Phenol index Total Kjeldahl Nitrogen Total Phosphorous	Acetic, Citric, Gluconic, Lactic, Malic and Tartric Acid Alpha Amino Nitrogen Ammonia Anthocyanin Arginine Copper Fructose Free SO2 Glucose Glycerol Iron Polyphenols Total Acidity Total SO2 Urea

Other parameters on demand. Kits supplied in liquid format in vials adapted to your Smartchem. 1 to 3 reagents necessary according the parameter.

THE AMS ALLIANCE DISCRETE ANALYZERS RANGE



Smartchem 140 Up to 140 tests/hour



Smartchem 200 Up to 200 tests/hour



SMARTCHEM

Smartchem 300 Up to 300 tests/hour



TECHNICAL SPECIFICATIONS

DESCRIPTION	Fully automated, direct read, floor sta	and, random access, Discrete analyzer			
ASSAY TYPE	End Point, Kinetic, Differential, Bi chromatic assays up to 999 tests, fully programmable				
ON-BOARD REAGENTS	Up to 100 reagents can be loaded into the refrigerated compartment and distributed on: - 5 removable racks (8 main positions) + 5 satellites - 100/50/25/5 ml containers - Reagents barcode reading (option)				
Samples loading	 Up to 200 samples on board 2 independent carousels hosting 20 racks of 10 positions each, for continuous sample loading Cups and/or tubes of 10-13 mm (max height 100 mm) Sample tube caps for sensitive application Samples barcode reading (option) 				
URGENT SAMPLE ANALYSIS	allowed at any time				
THROUGHPUT	Up to 600 tests/hour including ISE (pH, conductivity, redox)				
READING SYSTEM	Direct photometry				
MINIMUM REACTION VOLUME	180 μl				
	Photometer	direct reading with Reference channel			
	Wavelength	12 available wavelengths from 340 nm to 880 nm			
OPTICS SYSTEM	Optical Path	6 or 10 mm			
	Linear range	0.001 - 4.200 Abs			
Sampling system	 3 sampling arms 3 Shock sensors 1 blockage sensor 3 Capacitive liquid level sensors 2 positions for stirrer/mixer Automatic probe washing Sample rerun with or without dilution Limited water consumption (less than 7 liters/hour) 				
DILUTOR	Automatic dilution in pre and post-dilution				
REACTION PLATE	 85 cuvettes divided into 5 racks Automatic cuvette washing Continuous monitoring of cuvette quality Incubation temperature control at 37°C +/- 0,3 Separate waste lines to minimize hazard disposable management Use of washable cuvettes (Green technology) 				
DIMENSIONS	96 x 70 x 120 cm (W x D x H)				
WEIGHT	115 Kg				
OPTIONAL ACCESSORIES	- ISE module for pH, conductivity & Redox - Nitrate module with new reductor - Barcode reader for reagents - Barcode reader for samples				
Complementary products	Digestion system: Easy digest Dis	stillation system: Easy Dist Ready to use reagents			
AMS Alliance www.amsalliance.com	ORDER ONLI reagents and cons http://webshop.amsal				